



Main sol. 738

NH ₄ Cl	1.00	g
Na ₂ HPO ₄ x 2 H ₂ O	0.42	g
Na ₂ HPO ₄ x H ₂ O	0.18	g
MgCl ₂ x 6 H ₂ O	0.10	g
CaCl ₂ x 2 H ₂ O	0.10	g
KCl	0.10	g
FeSO ₄ x 7 H ₂ O	0.04	g
Wolin's vitamin solution	10.00	ml
Modified Wolin's mineral solution	10.00	ml
Distilled water	1000.00	ml
Cresol or Phenol	1.00	mM

1. Dissolve phosphates separately, Adjust pH to 6.5. Add cresol/phenol and vitamin solution after autoclaving.

2. (e.g.: 20 ml/l of a 4.7 g/l phenol-solution)