



Main sol. 172

Yeast extract (Difco)	1.00	g
Tryptone (Difco)	1.00	g
NaCl	24.70	g
KCl	0.70	g
MgSO ₄ x 7 H ₂ O	6.30	g
MgCl ₂ x 6 H ₂ O	4.60	g
CaCl ₂ x 2 H ₂ O	1.20	g
NaHCO ₃	0.20	g
Agar (Difco)	15.00	g
Distilled water	1000.00	ml

Adjust pH to 7.2. Sodium bicarbonate and calcium chloride are autoclaved separately, each in a small volume of distilled water.