



## **Pfennig's heterotrophic salts solution** (from medium 28)

Ammonium chloride	0.35	g
Ammonium acetate	0.25	g
Pyruvic acid sodium salt	0.25	g
Dextrose	0.25	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	0.50	g
KCl	0.35	g
KH <sub>2</sub> PO <sub>4</sub>	0.35	g
<b>Trace element solution SL-12 B</b>	1.00	ml
Distilled water	25.00	ml

Filter sterilize into sterile, gas-tight, 100 ml screw-cap bottle.