



Main sol. 1727

KH ₂ PO ₄ (0.41 g/l stock solution)	100.00	ml
K ₂ HPO ₄ (4.35 g/l stock solution)	3.90	ml
CaCl ₂ x 2 H ₂ O (3.6 g/l stock solution)	3.60	ml
MgSO ₄ x 7 H ₂ O (7.5 g/l stock solution)	3.33	ml
Fe-EDTA complex	10.00	ml
Micronutrient solution (after H. Gaffron)	0.06	ml

1. Make up to 1000 ml with MilliQ water.
2. Autoclave at 121 °C for 20 min.