



Trace metal solution

EDTA	50.00	g
ZnSO ₄ x 7 H ₂ O	22.00	g
CaCl ₂	5.54	g
MnCl ₂ x 4 H ₂ O	5.06	g
FeSO ₄ x 7 H ₂ O	4.99	g
(NH ₄) ₆ Mo ₇ O ₂₄ x 4 H ₂ O	1.10	g
CuSO ₄ x 5 H ₂ O	1.57	g
CoCl ₂ x 6 H ₂ O	1.61	g
Distilled water	1000.00	ml

Adjust pH to 6.0 with KOH.