



1203: FASTIDIOUS ANAEROBE AGAR

This recipe contains strain-specific modifications for *Vermiculatibacterium agrestimuris* DSM 112226 *

Final pH: 7.0 - 7.4

Final volume: 1000 ml

Fastidious Anaerobe Agar (acumedia)	45.70	g
Deionized water	1000.00	ml
Horse blood, adjust if required	10.00	%

Allow to soak for 10 minutes, swirl to mix then sterilize by autoclaving at 121°C for 15 minutes. Cool to 47°C then aseptically add 5-10% of sterile defibrinated horse blood, mix well and pour plates. pH: 7.2 ± 0.2

* solid medium, add. 5% horse blood