

## 1393: HALOMONAS WANG MEDIUM

This recipe contains strain-specific modifications for *Marinobacter sp.* DSM 22447 \*

Final pH: \* 8.5

Final volume: 1000 ml

### **Mineral Salts Solution**

Casamino acids	5.00	g
Yeast extract	5.00	g
NaCl	50.00	g
Distilled water	1000.00	ml
Agar, for solid medium	20.00	g

Adjust to pH 9.0 - 9.5 and sterilize by autoclaving. The medium may be solidified by adding 20.0 g/l agar.

\* Adjust pH of final medium to 8.5.