

195c: DESULFOBACTER SP. MEDIUM (LACTATE)

This recipe contains strain-specific modifications for *Paucidesulfovibrio capillatus* DSM 14982 *

Final pH: 7.1 - 7.4

Final volume: 1003 ml

Solution A	952.00	ml
Solution B	30.00	ml
Solution C	10.00	ml
Solution D	1.00	ml
Solution E	10.00	ml

1. Solution A is sparged with 80% N₂ and 20% CO₂ gas mixture to reach a pH below 6 (at least 30 min), then distributed under the same gas atmosphere in anoxic Hungate-type tubes or serum vials and autoclaved. Solutions C and E are autoclaved separately under 100% N₂ gas. Solution B is autoclaved under 80% N₂ and 20% CO₂ gas atmosphere. Solution D is prepared under 100% N₂ gas atmosphere and sterilized by filtration. To complete the medium appropriate amounts of solutions B to E are added to the sterile solution A in the sequence as indicated. Final pH of the medium should be 7.1 - 7.4.

2. Note: Addition of 10 - 20 mg sodium dithionite per liter (e.g. from 5% (w/v) solution, freshly prepared under N₂ and filter-sterilized) may stimulate growth at the beginning. For transfers use 5 - 10% (v/v) inoculum.

* Supplement medium with 5.00 g/l Na₂S₂O₃ x 5 H₂O added to the autoclaved medium from an anoxic stock solution sterilized by filtration.

Solution A

Na ₂ SO ₄	3.00	g
KH ₂ PO ₄	0.20	g
NH ₄ Cl	0.30	g
NaCl	21.00	g
MgCl ₂ x 6 H ₂ O	3.00	g
KCl	0.50	g
CaCl ₂ x 2 H ₂ O	0.15	g
Trace element solution SL-10	1.00	ml
Selenite-tungstate solution	1.00	ml
Sodium resazurin (0.1% w/v)	0.50	ml
Na ₂ S ₂ O ₃ x 5 H ₂ O	5.00	g
Distilled water	950.00	ml

Solution B

Na ₂ CO ₃	1.50	g
---------------------------------	------	---

195c: DESULFOBACTER SP. MEDIUM (LACTATE)

Distilled water	30.00	ml
-----------------	-------	----

Solution C

Na-L-lactate	2.50	g
Distilled water	10.00	ml

Solution D

Wolin's vitamin solution (10x)	1.00	ml
---------------------------------------	------	----

Solution E

Na ₂ S x 9 H ₂ O	0.40	g
Distilled water	10.00	ml

Selenite-tungstate solution (from medium 385)

NaOH	0.50	g
Na ₂ SeO ₃ x 5 H ₂ O	3.00	mg
Na ₂ WO ₄ x 2 H ₂ O	4.00	mg
Distilled water	1000.00	ml

Trace element solution SL-10 (from medium 320)

HCl (25%)	10.00	ml
FeCl ₂ x 4 H ₂ O	1.50	g
ZnCl ₂	70.00	mg
MnCl ₂ x 4 H ₂ O	100.00	mg
H ₃ BO ₃	6.00	mg
CoCl ₂ x 6 H ₂ O	190.00	mg
CuCl ₂ x 2 H ₂ O	2.00	mg
NiCl ₂ x 6 H ₂ O	24.00	mg
Na ₂ MoO ₄ x 2 H ₂ O	36.00	mg
Distilled water	990.00	ml

First dissolve FeCl₂ in the HCl, then dilute in water, add and dissolve the other salts. Finally make up to 1000.00 ml.

Wolin's vitamin solution (10x) (from medium 120)

Biotin	20.00	mg
Folic acid	20.00	mg
Pyridoxine hydrochloride	100.00	mg
Thiamine HCl	50.00	mg
Riboflavin	50.00	mg



195c: DESULFOBACTER SP. MEDIUM (LACTATE)

Nicotinic acid	50.00	mg
Calcium D-(+)-pantothenate	50.00	mg
Vitamin B ₁₂	1.00	mg
p-Aminobenzoic acid	50.00	mg
(DL)-alpha-Lipoic acid	50.00	mg
Distilled water	1000.00	ml