

31: ALKALINE NUTRIENT AGAR

This recipe contains strain-specific modifications for *Alkalicoccobacillus gibsonii* DSM 9728 *

Final pH: 9.7

Final volume: 1000 ml

Peptone	5.00	g
Meat extract	3.00	g
Agar, if necessary	15.00	g
Distilled water	900.00	ml

1. After sterilization add sterile 1 M Na-sesquicarbonate solution (1 ml in 10 ml) to achieve a pH of 9.7

Na-sesquicarbonate solution	100.00	ml
NaCl	5.00	g/l

2. For Bacillus strains the addition of 10.0 mg $\text{MnSO}_4 \times \text{H}_2\text{O}$ is recommended for sporulation.

* With 0.5% NaCl

Na-sesquicarbonate solution (from medium 31)

NaHCO_3	4.20	g
Na_2CO_3 anhydrous	5.30	g
Distilled water	100.00	ml