

545: TRYPTONE SOYA BROTH (TSB)

This recipe contains strain-specific modifications for *Thalassospira* sp. DSM 17669 *

Final pH: 7.3

Final volume: 1000 ml

Casein peptone	17.00	g
Soy peptone	3.00	g
D(+)-Glucose	2.50	g
NaCl	5.00	g
K ₂ HPO ₄	2.50	g
NaCl	15.00	g/l
Distilled water	1000.00	ml

1. Adjust pH to 7.3.
2. The same as Caso Bouillon (Merck or Oxoid, Germany).

* Add 1.5% NaCl