

545: TRYPTONE SOYA BROTH (TSB)

This recipe contains strain-specific modifications for *Aquamicrobium aerolatum* DSM 21857 *

Final pH: * 7.2

Final volume: 1000 ml

Casein peptone	17.00	g
Soy peptone	3.00	g
D(+)-Glucose	2.50	g
NaCl	5.00	g
K ₂ HPO ₄	2.50	g
Distilled water	1000.00	ml

1. Adjust pH to 7.3.
2. The same as Caso Bouillon (Merck or Oxoid, Germany).

* PH 7.2