

549: CM + YE MEDIUM

This recipe contains strain-specific modifications for *Virgibacillus campisalis* DSM 25790 *

Final volume: 1000 ml

Casamino acids (vitamin free, Difco)	7.50	g
Yeast extract	10.00	g
MgSO ₄ x 7 H ₂ O	20.00	g
Na ₃ -citrate	3.00	g
KCl	2.00	g
NaCl	50.00	g/L
Fe solution	1.00	ml
Agar	15.00	g
Distilled water	1000.00	ml

pH 7.4 (NaOH)

* With 50g/L NaCl only

Fe solution

FeSO ₄ x 7 H ₂ O	10.00	g
HCl (0.001 M)	1000.00	ml

Prepare 4.98% FeSO₄ x 7 H₂O (= 10 g Fe in 1l 0.001 M HCl)