



## 550: CYC-MEDIUM (modified following Cross and Attwell,1973)

This recipe contains strain-specific modifications for *Tuberibacillus calidus* DSM 17572 \*

Final pH: 7.2

Final volume: 1000 ml

Czapek Dox agar (Merck)	48.00	g
Yeast extract (Oxoid)	2.00	g
Casamino acids (Difco)	6.10	g
Tryptophan	0.02	g
Distilled water	1000.00	ml
Agar	20.00	g
<b>MnSO<sub>4</sub></b>	<b>10.00</b>	<b>mg/l</b>

pH 7.2

\* Plus MnSO<sub>4</sub> (10µg/ml)