

550: CYC-MEDIUM (modified following Cross and Attwell,1973)

This recipe contains strain-specific modifications for *Thermoactinomyces vulgaris* DSM 43368 *

Final pH: * 8.0

Final volume: 1000 ml

Czapek Dox agar (Merck)	48.00	g
Yeast extract (Oxoid)	2.00	g
Casamino acids (Difco)	6.10	g
Tryptophan	0.02	g
Distilled water	1000.00	ml
Agar	20.00	g

pH 7.2

* pH 8