

618: TWEEN 80A MEDIUM

Final pH: 7.2

Final volume: 1000 ml

Casitone (Difco)	5.00	g
Yeast extract (Difco)	3.00	g
MgSO ₄ x 7 H ₂ O	2.00	g
Distilled water	1000.00	ml
Agar, for solid medium	15.00	g/l
Tween 80 (sterile)	1.00	%

For solid media add 15.0 g agar. Dissolve ingredients and adjust the pH to 7.2. Sterilize by autoclaving. After sterilization add sterile Tween 80 to a final concentration of 1.0%.