

## 65: GYM STREPTOMYCES MEDIUM

This recipe contains strain-specific modifications for *Nocardiopsis dassonvillei* DSM 44686 \*

Final pH: \* 7.2

Final volume: 1000 ml

|                                 |         |    |
|---------------------------------|---------|----|
| Glucose                         | 4.00    | g  |
| Yeast extract                   | 4.00    | g  |
| Malt extract                    | 10.00   | g  |
| CaCO <sub>3</sub> , if required | 2.00    | g  |
| Agar                            | 20.00   | g  |
| Distilled water                 | 1000.00 | ml |

Adjust pH to 7.2 before adding agar. Delete CaCO<sub>3</sub> if liquid medium is used.

\* Adjust pH of final medium to 7.2.