## **Microorganisms**



## 830: R2A MEDIUM

This recipe contains strain-specific modifications for Ramlibacter tataouinensis DSM 14655 \*

Final pH: 7.2 Final volume: 1000 ml

Yeast extract	0.50	g
Proteose peptone (Difco no. 3)	0.50	g
Casamino acids	0.50	g
Glucose	0.50	g
Starch (soluble)	0.50	g
Na-pyruvate	0.30	g
K <sub>2</sub> HPO <sub>4</sub>	0.30	g
$MgSO_4 \times 7 H_2O$	0.05	g
Agar	15.00	g
Distilled water	1000.00	ml

Final pH 7.2; adjust with crystalline  $K_2$ HPO<sub>4</sub> or KH<sub>2</sub>PO<sub>4</sub> before adding agar. Add agar, heat medium to boiling to dissolve agar, and autoclave for 15 min at 121°C.

\* Rehydrate and grow lyophilized cells from the ampoule in 5 ml liquid broth 1. Subsequent subculturing may be carried out in liquid medium or on agar plates