

## 92: TRYPTICASE SOY YEAST EXTRACT MEDIUM

This recipe contains strain-specific modifications for *Microlunatus sp.* DSM 21973 \*

Final pH: \* 6.0 - 7.7

Final volume: 1000 ml

Trypticase soy broth	30.00	g
Yeast extract	3.00	g
Agar, if necessary	15.00	g
Distilled water	1000.00	ml

Adjust pH to 7.0 - 7.2.

\* Adjust pH of final medium to 6.0 - 7.7.