

381: LB (Luria-Bertani) MEDIUM

This recipe contains strain-specific modifications for *Chromobacterium vaccinii* DSM 25150 *

Final pH: 7.0

Final volume: 1000 ml

Tryptone	10.00	g
Yeast extract	5.00	g
NaCl	10.00	g
Agar	20.00	g
Distilled water	1000.00	ml

Adjust pH to 7.0.

* Rehydrate and grow lyophilized cells in liquid broth 381. Subsequent subculturing may be carried out in liquid medium or on agar plates