



545: TRYPTONE SOYA BROTH (TSB)

This recipe contains strain-specific modifications for *Francisella orientalis* DSM 21254 *

Final pH: 7.3

Final volume: 1000 ml

Casein peptone	17.00	g
Soy peptone	3.00	g
D(+)-Glucose	2.50	g
NaCl	5.00	g
K ₂ HPO ₄	2.50	g
Distilled water	1000.00	ml

1. Adjust pH to 7.3.

2. Note: The following brands are used at the DSMZ: Difco™ Trypton-Soja-Bouillon (Soja-Caseinpepton-Bouillon) BD 211825. BBL™ Trypticase Soja-Bouillon (Caseinpepton-Soja-Bouillon) BD 211768. Oxoid Tryptone Soy Broth (Casein soya bean digest medium) CM0129 Merck Tryptic Soy Broth (CASO-Bouillon) 105459

* plus 5g/L Glucose plus 1g/L L-Cystein